



## AVAILABLE COMBO COOLERS & ESTIMATED DOSAGE CAPACITY

The following are cooler-in-a-cooler combinations designed for double-insulating swine semen shipments, from smallest to largest:

INSIDE COOLER (wall thickness)	OUTSIDE COOLER (wall thickness)	ESTIMATED DOSES	
		BAGS	BOTTLES/TUBES
SK2 (1.375")	TK9 (1.5")	2-4	2
SK3 (1.375")	TK9 (1.5")	2-6	2-4
S6x9 (.75")	TK9 (1.5")	8-10	4-6
SK5 (1.375")	SL810 (1.5")	10-12	8-10
TK4 (1.5")	SL10 (1.5")	20-25	15-20
TK4 (1.5")	PC18 (1.875")	20-25	15-20
SG4.5 (1.5")	PC18 (1.875")	30-35	20-25
SG4.5 (1.5")	SL916 (1.5")	30-35	20-25
TK6 (1.5")	TK20 (1.5")	35-40	25-30
TK6 (1.5")	SL10 (1.5")	35-40	25-30
TK7 (1")	TK20 (1.5")	40-45	30-35
TK8 (1.5")	PC30 (1.875")	45-50	35-40
TK8 (1.5")	SL1116 (1.5")	45-50	35-40
TK9 (1.5")	XC34 (1.5")	60-70	50-60
SL610 (1.5")	SL922 (1.5")	70-80	55-65
SL10 (1.5")	SL22 (1.5")	110-125	80-100
LC8 (1.5")	PC53 (1.5")	150-180	130-140
LC8 (1.5")	PC64 (1.5")	150-180	130-140
LC8 (1.5")	PC75 (1.5")	150-180	130-140
SL922 (1.5")	PC90 (1.5")	200-225	150-175
PC850 (1.5")	PC80 (1.5")	250	200
PC1050 (1.5")	PC90 (1.5")	300	250

The estimated dose capacity quantities are based on feedback from our current customers, and may vary for your application. Samples are available for test packing to insure proper dosage fit. Visit the *Insulated Shipping Cooler* page of our website for specifications of the coolers listed above.